

Human BMF protein (Recombinant) (T7-Tag) (STJP016470)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing, 1mM DTT, 10% Glycerol
Molecular Weight	15.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	90427
Gene Symbol	BMF
UniProt ID	BMF_HUMAN
Immunogen Region	1-129aa
Immunogen Sequence	MASMTGGQQM GRGSHMEPSQ CVEELEDVDF QPEDGEPVTQ PGSLLSADLF AQSLLDLCPLS RLQLFPLTHC CGPGLRPTSQ EDKATQTLSP ASPSQGVMLP CGVTEEPQRL FYAPAEPKSC VVADPPLPAQ PCFEWRREQE RGRP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human BMF protein 3µg under reducing conditions Coomassie blue stain