

Human SAP18 protein (Recombinant) (His-Tag) (STJP016483)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 30% Glycerol, 1mM DTT
Molecular Weight	19.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10284
Gene Symbol	SAP18
UniProt ID	SAP18_HUMAN
Immunogen Region	20-172aa
Immunogen Sequence	MAVESRVLTQE EIKKEPEKPI DREKTCPLLL RVFTTNNGRH HRMDEFSGN VPSELQIYT WMDATLKELT SLVKEVYPEA RKKGTHFNFA IVFTDVKRPG YRVKEIGSTM SGRKGTDDSM TLQSQKFQIG DYLDIAITPP NRAPPPSGRM RPY

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human SAP18 protein under reducing conditions stained with Coomassie blue