

Human CRABP1 protein (Recombinant) (No-Tag) (STJP016494)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	15.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1381
Gene Symbol	CRABP1
UniProt ID	RABP1_HUMAN
Immunogen Region	1-137aa
Immunogen Sequence	MPNFAGTWKM RSSENFDELL KALGVNAMLK KVAVAAASKP HVEIRQDGDQ FYIKTSTTVR TTEINFKVGE GFEETVDGR KCRSLATWEN ENKIHCTQL LEGDGPITYW TRELANDELI LTFGADDVVC TRIYVRE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human CRABP1 protein under reducing conditions stained with Coomassie blue