

Human RAC2 protein (Recombinant) (His-Tag) (STJP016507)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol 2mM EDTA, 2mM DTT, 0.1 M NaCl.
Molecular Weight	23.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5880
Gene Symbol	RAC2
UniProt ID	RAC2_HUMAN
Immunogen Region	1-189aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MQAIKCVVVG DGAVGKTCLL ISYTTNAFPG EYIPTVFDNY SANVMVDSKP VNLGLWDTAG QEDYDRLRPL SYPQTDVFLI CFSLVSPASY ENVRAKWFPE VRHHCPTPI ILVGTKLDLR DDKDTIEKLEK EKKLAPITYP QGLALAKEID SVKYLECSAL TORGLKTVFD EAIRAVLCPQ PTRQQKRC

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RAC2 protein under reducing conditions stained with Coomassie blue