

## Human CCL5/RANTES protein (Recombinant) (No-Tag) (STJP016516)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 10mM Sodium Citrate buffer (pH 3.5) containing 10% Glycerol
Molecular Weight	7.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6352</a>
Gene Symbol	<a href="#">CCL5</a>
UniProt ID	<a href="#">CCL5_HUMAN</a>
Immunogen Region	24-91aa
Immunogen Sequence	MSPYSSDTP CCFAYIARPL PRAHIKEYFY TSGKCSNPAV VVTRKNRQV CANPEKKWVR EYINSLEMS

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human CCL5 RANTES protein under reducing conditions