

## Human CX3CL1/Fractalkine protein (Recombinant) (His-Tag) (STJP016522)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) 10% Glycerol
Molecular Weight	10.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6376</a>
Gene Symbol	<a href="#">CX3CL1</a>
UniProt ID	<a href="#">X3CL1_HUMAN</a>
Immunogen Region	25-100aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MQHHGVTKCN ITC SKMTSKI PVALLIHYQQ NQASCGKRAI ILETRQHRLF CADPKEQWVK DAMQHLD RQA AALTRNG

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human CX3CL1 Fractalkine protein under reducing conditions  
 Coomassie blue stain