

Human Galectin-1/LGALS1 protein (Recombinant) (No-Tag) (STJP016527)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 10% Glycerol
Molecular Weight	14.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3956
Gene Symbol	LGALS1
UniProt ID	LEG1_HUMAN
Immunogen Region	1-135aa
Immunogen Sequence	MACGLVASNL NLKPGELRV RGEVAPDAKS FVLNLGKDSN NLCLHFNPRF NAHGDANTIV CNSKDGGAWG TEQREAVFPF QPGSVAEVC I TFDQANLTVK LPDGYEFKFP NRLNLEAINY MAADGDFKIK CVAFD

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Galectin-1 LGALS1 protein under reducing conditions stained with Coomassie blue