

Human CDC42 protein (Recombinant) (T7-Tag) (STJP016537)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% Glycerol, 2mM EDTA
Molecular Weight	22.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	998
Gene Symbol	CDC42
UniProt ID	CDC42_HUMAN
Immunogen Region	1-188aa
Immunogen Sequence	MQTIK CVVVG DGAVG KTCLLSYTT NKFPSEYVPT VFDNYAVTVM IGGEPYTLGL FDTAGQEDYD RLRPLSYPT DVFLVCFSVV SPSSFENVKE KVVPEITHHC PKTPFLLVGT QIDLRDDPST IEKLAKNKQK PITPETAELK ARDLKAVKYV ECSALTQKGL KNVFDEAILA ALEPPEPKKS RRC

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human CDC42 protein under reducing conditions stained with Coomassie blue