

Human Aldehyde Dehydrogenase 1-A1/ALDH1A1 protein (Recombinant) (No-Tag) (STJP016547)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 50mM Tris-HCl buffer (pH 7.5) containing 10% Glycerol
Molecular Weight	54.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	216
Gene Symbol	ALDH1A1
UniProt ID	AL1A1_HUMAN
Immunogen Region	1-501aa
Immunogen Sequence	MSSSGTDPDLP VLLTDLKIY TKIFINNEWH DSVSGKKFPV FNPATEEELC QVEEGDKEDV DKAVKAARQA FQIGSPWRTM DASERGRLLY KLADLIERDR LLLATMESMN GKKLYSNAYL NDLAGCIKTL RYCAGWADKI QGR TIPIDGN FFTYTRHEPI GVCQIIPWN FPLVMLIWKI GPALSCGNTV VVKPAEQTPL TALHVASLIK EAGFPPGVVN IVPGYGPT

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human ALDH1A1 protein under reducing conditions Coomassie blue stain