

Human BAG3 protein (Recombinant) (His-Tag) (STJP016582)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM EDTA, 0.1mM PMSF, 10% Glycerol
Molecular Weight	63.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9531
Gene Symbol	BAG3
UniProt ID	BAG3_HUMAN
Immunogen Region	1-575aa
Immunogen Sequence	MSAATHSPMM QVASNGDRD PLPPGWEIKI DPQTGWPFV DHNSRTTWN DPRVPSEGPK ETPSSANGPS REGSRLPPAR EGHVPVQLR PGYIPIVLH EGAENRQVHP FHVYPQGMQ RFRTEAAAAA PQRSQSPLRG MPETTQPKQ CGQVAAAAA QPPASHGPER SQSPAASDCS SSSSSASLPS SGRSSLGSHQ LPRGYI

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human BAG3 protein under reducing conditions Coomassie blue stain