

Human PTS protein (Recombinant) (His-Tag) (STJP016594)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% Glycerol
Molecular Weight	18.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5805
Gene Symbol	PTS
UniProt ID	PTPS_HUMAN
Immunogen Region	1-145aa
Immunogen Sequence	MSTEGGGRRRC QAQVSRRI SF SASHRLYSKF LSDEENLKLF GKCNNPNHG HNYKVVVTVH GEIDPATGMV MNLADLKKYM EEAIMQPLDH KNLDMDVPYF ADVVSTTEN VAVYIWDNLQK VLPVGVLYKV KYYETDNNIV VYKGE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human PTS protein under reducing conditions with Coomassie blue stain