

Human Annexin A4/ANXA4 protein (Recombinant) (No-Tag) (STJP016605)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol 0.2M NaCl
Molecular Weight	36 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	307
Gene Symbol	ANXA4
UniProt ID	ANXA4_HUMAN
Immunogen Region	1-321aa
Immunogen Sequence	MAMATKGGTV KAASGFNAME DAQTLRKAMK GLGTDEDAII SVLAYRNTAQ RQEIRTAYKS TIGRDLIDDL KSELSGNFEQ VIVGMMTPPTV LYDVQELRRA MKGAGTDEGC LIEILASRTP EEIRRIISQTY QQQYGRSLED DIRSDTSFMF QRVLVLSAG GRDEGNYLDD ALVRQDAQDL YEAGEKKWGT DEVKFLTVLC SRNRNHLHVF DEYKRISQK DIEQSIKS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE of recombinant human Annexin A4 ANXA4 protein under reducing conditions