

Human Annexin A2/ANXA2 protein (Recombinant) (His-Tag) (STJP016609)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol 0.1M NaCl
Molecular Weight	42.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	302
Gene Symbol	ANXA2
UniProt ID	ANXA2_HUMAN
Immunogen Region	1-339aa
Immunogen Sequence	MST VHEILCKLSL EGDHSTPPSA YGSVKAYTNF DAERDALNIE TAIKTKGVDE VTIVNILTNR SNAQRQDI AFAYQRRTKKEL ASALKSALSG HLETVILGLL KTPAQYDASE LKASMKGLGT DEDSLIEIIC SRTNQELQEI NRVYKEMYKT DLEKDIISDT SGDFRKL MVA LAKGRR AEDG SVIDYE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human Annexin A2 ANXA2 protein under reducing conditions