

Human CDK2AP1 protein (Recombinant) (His-Tag) (STJP016621)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 50% Glycerol, 2mM DTT
Molecular Weight	16.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8099
Gene Symbol	CDK2AP1
UniProt ID	CDKA1_HUMAN
Immunogen Region	1-115aa
Immunogen Sequence	MRGSHHHHHH GMASMTGGQQ MGRDLYDDDD KDRWGSMSY KPNLAHMPA AALNAAGSVH SPSTSMATSS QYRQLSDYG PPSLGYTQGT GNSQVPQSKY AELLAIEEL GKEIRPTYAG SKSAMERLKR GIIHARGLVR ECLAETERNA RS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human CDK2AP1 protein under reducing conditions stained with Coomassie blue