

Human IVD protein (Recombinant) (His-Tag) (STJP016628)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% Glycerol
Molecular Weight	45.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3712
Gene Symbol	IVD
UniProt ID	IVD_HUMAN
Immunogen Region	33-426aa
Immunogen Sequence	MHSLLPVDDA INGLSEEQRQ LRQTMKFLQ EHLAPKAEI DRSNEFKNLR EFWKQLGNLG VLGITAPVQY GGSGGLGYLEH VLMEEISRA SGAVGLSYGA HSNLCINQLV RINGNEAQKEK YLPKLISGEY IGALAMSEPN AGSDVSMKL KAEKKGNYI LNGNKFWITN GPDADVLIVY AKTDLAAPVA SRGITAFIVE KGMPGF

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human IVD protein 3µg under reducing conditions Coomassie blue stain