

Human FKBP13/FKBP2 protein (Recombinant) (No-Tag) (STJP016647)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT
Molecular Weight	13.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2286
Gene Symbol	FKBP2
UniProt ID	FKBP2_HUMAN
Immunogen Region	22-142aa
Immunogen Sequence	MATGAEGKRRK LQIGVKRVD HCPKSRKGD VLHMHYTGKL EDGTEFDSSL PQNQPFVFSL GTGQVIKQWD QLLGMCEGE KRKLVIPSEL GYGERGAPPK IPGGATLVFE VELLKIERRT EL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel of recombinant human FKBP13/FKBP2 protein under reducing conditions
 Coomassie blue stain