

## Human BAG1 protein (Recombinant) (No-Tag) (STJP016649)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | SDS-PAGE |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                     |
|---------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 1 mg/mL                                                                                                                                             |
| Formulation         | Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 100mM NaCl, 10% Glycerol                                                                         |
| Molecular Weight    | 25.9 kDa                                                                                                                                            |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                                               |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">573</a>                                                                                                                                                                                                                                           |
| Gene Symbol        | <a href="#">BAG1</a>                                                                                                                                                                                                                                          |
| UniProt ID         | <a href="#">BAG1_HUMAN</a>                                                                                                                                                                                                                                    |
| Immunogen Region   | 1-230aa                                                                                                                                                                                                                                                       |
| Immunogen Sequence | MNRSQEVTRD EESTRSEEV T REEMAAAGLT VTVTHSNEKH DLHVTSQQGS SEPVVQDLAQ VVEEVIGVPQ SFQKLIFKGK SLKEMETPLS ALGIQDGCVRV MLIGKKNSPQ EEVELKKLKH LEKSVEKIAD QLEELNKELT GIQQGF LPKD LQAEALCKLD RRVKATIEQF MKILEEIDTL ILPENFKDSR LKRKGLVKKV QAFLAECDTV EQNICQETER LQSTNFAL |

### ADDITIONAL INFORMATION

|      |                                                                                                                         |
|------|-------------------------------------------------------------------------------------------------------------------------|
| Note | STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS. |
|------|-------------------------------------------------------------------------------------------------------------------------|



SDS-PAGE gel of recombinant human BAG1 protein under reducing conditions stained with Coomassie blue