

Human UBE2T protein (Recombinant) (His-Tag) (STJP016669)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 1mM DTT
Molecular Weight	23.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	29089
Gene Symbol	UBE2T
UniProt ID	UBE2T_HUMAN
Immunogen Region	1-197aa
Immunogen Sequence	MQRASRLKRE LHMLATEPPP GITCWQDKDQ MDDLRAQILG GANTPYEKG V FKLEVIIPER YPFEPPIQIRF LTPIYHPNID SAGRICLDVL KLPPKGAWRP SLNIATVLTSL IQLLMSEPNP DDPLMADISS EFKYNKPAFL KNARQWTEKH ARQKQKADEE EMLDNLPEAG DSRVHNSTQK RKASQLVGIE KKFHPDVLEH HHHHHH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human UBE2T protein under reducing conditions stained with Coomassie blue