

Human SGTA protein (Recombinant) (His-Tag) (STJP016673)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 1mM DTT
Molecular Weight	35.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6449
Gene Symbol	SGTA
UniProt ID	SGTA_HUMAN
Immunogen Region	1-313aa
Immunogen Sequence	MDNKKRLAYA IIQFLHDQLR HGGSSDAQE SLEVAIQCLE TAFGVTVEDS DLALPQTLPE IFEEAAATGKE MPQDLRSPAR TPPSEEDSAE AERLKTEGNE QMKVENFEAA VHFYKGKIEL NPANAVYFCN RAAAYSGLGN YAGAVQDCER AICIDPAYSK AYGRMGLALS SLNKHVEAVA YYKKALELDP DNETYKSNLK IAEKLRAP SPTGGVGSFD IAGLLNNP

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human SGTA protein 3µg under reducing conditions stained with Coomassie blue