

Human CUTA protein (Recombinant) (His-Tag) (STJP016681)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% Glycerol
Molecular Weight	17.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51596
Gene Symbol	CUTA
UniProt ID	CUTA_HUMAN
Immunogen Region	33-179aa
Immunogen Sequence	MRLLLLPRVL LTMASGSPPT QPSPASDSGS GYVPGSVSAA FVTCPNKVA KEIARAVVEK RLAACVNLIP QITSIYEWKG KIEEDSEVLM MIKTQSSLVP ALTDVFRSVH PYEVAEVIAL PVEQGNFPYL QWVRQVTESV SDSITVLP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human CUTA protein 3µg under reducing conditions stained with Coomassie blue