


E.coli trxA protein (Recombinant) (His-Tag) (STJP016689)

GENERAL INFORMATION		
	Product Type	Proteins
	Applications	SDS-PAGE/Enzyme Activity
	Host / Source	E.coli
PRODUCT PROPERTIES		
	Concentration	0.5 mg/mL
	Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT
	Molecular Weight	12.8 kDa
	Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION		
	Immunogen Region	2-109aa
	Immunogen Sequence	MHHHHHHMGS DKIIHLTDDS FDTDVLKADG AILVDFWAEW CGPCKMIAP I LDEIADEYQG KLTVAKLNID QNPGTAPKYG IRGIPTLLLF KNGEVAATKV GALSKGQLKE FLDANLA
ADDITIONAL INFORMATION		
	Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant E.coli trxA protein 3µg under reducing conditions Coomassie blue stain