

## Human RAB27A protein (Recombinant) (His-Tag) (STJP016696)

### GENERAL INFORMATION

|               |          |
|---------------|----------|
| Product Type  | Proteins |
| Applications  | SDS-PAGE |
| Host / Source | E.coli   |

### PRODUCT PROPERTIES


|                     |                                                                                                                                                    |
|---------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| Concentration       | 1 mg/mL                                                                                                                                            |
| Formulation         | Liquid in phosphate-Buffered Saline (pH 7.4) containing 20% Glycerol                                                                               |
| Molecular Weight    | 27 kDa                                                                                                                                             |
| Storage Instruction | For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles. |

### TARGET INFORMATION

|                    |                                                                                                                                                                                                                                  |
|--------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Gene ID            | <a href="#">5873</a>                                                                                                                                                                                                             |
| Gene Symbol        | <a href="#">RAB27A</a>                                                                                                                                                                                                           |
| UniProt ID         | <a href="#">RB27A_HUMAN</a>                                                                                                                                                                                                      |
| Immunogen Region   | 1-221aa                                                                                                                                                                                                                          |
| Immunogen Sequence | MSDGDYDYLI KFLALGDSGV GKTSVLYQYT DGKFNSKFIT TVGIDFREKR VVYRASGPDG ATGRGQRIHL QLWDTAGQER FRSLTTAFFR DAMGFLLFD LTNEQSFLNV RNWISQLQMH AYCENPDIIVL CGNKSDLEDQ RVVKEEEAIA LAEKYGIPYF ETSANGTNI SQAIEMLLDL IMKRMECVD KSWIPEGVVR SNGHAS |

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



Recombinant human RAB27A protein SDS-PAGE gel under reducing conditions Coomassie blue stain