

Human RAB11A protein (Recombinant) (His-Tag) (STJP016697)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% Glycerol, 1mM EDTA
Molecular Weight	28.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8766
Gene Symbol	RAB11A
UniProt ID	RB11A_HUMAN
Immunogen Region	1-213aa
Immunogen Sequence	MGTR DDEYDYLFKV VLGDSGVGK SNLLSRFTRN EFNLESKSTI GVEFATRSIQ VDGKTIKAQI WDTAGQERYR AITSAYYRGA VGALLVYDIA KHLTYENVER WLKELRDHAD SNIVIMLVGN KSDLRHLRAV PTDEARAF AE KNGLSFIETS ALDSTNVEAA FQTILTEIYR IVSQKQMSDR RENDMS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RAB11A protein under reducing conditions Coomassie blue stain