

Human FKBP36/FKBP6 protein (Recombinant) (His-Tag) (STJP016700)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Enzyme Activity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 40% Glycerol, 1mM EDTA, 0.2M NaCl
Molecular Weight	39.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8468
Gene Symbol	FKBP6
UniProt ID	FKBP6_HUMAN
Immunogen Region	1-327aa
Immunogen Sequence	MGGSALNQGVLV LEGDDAPGQS LYERLSQRML DISGDRGVLK DVIREGAGDL VAPDASVLVK YSGYLEHMDR PFDSNYFRKT PRLMKLGEDI TLWGMELGLL SMRRGELARF LFKPNYAYGT LGCPPLIPPN TTVLFEIELL DFLDCAESDK FCALSAEQD QFPLQKVLKV AATEREFGNY LFRQNRFYDA KVRYSKRALLL LRRRSA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human FKBP36
FKBP6 protein 3µg under reducing conditions
Coomassie blue stain