

Human ACBP/DBI protein (Recombinant) (His-Tag) (STJP016707)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 10% Glycerol
Molecular Weight	12.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1622
Gene Symbol	DBI
UniProt ID	ACBP_HUMAN
Immunogen Region	1-87aa
Immunogen Sequence	MSQAEFEKAA EEVRHLKTKP SDEEMLFYIG HYKQATVGDI NTERPGMLDF TGKAKWDawn ELKGTSKEDA MKAYINKVEE LKKKYGI

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human ACBP DBI protein under reducing conditions stained with Coomassie blue