

Human Ubc2/UBE2B protein (Recombinant) (His-Tag) (STJP016709)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 0.1M NaCl
Molecular Weight	21.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7320
Gene Symbol	UBE2B
UniProt ID	UBE2B_HUMAN
Immunogen Region	1-152aa
Immunogen Sequence	MSTP ARRRLMRDFK RLQEDPPVGV SGAPSENNIM QWNAVIFGPE GTPFEDGTFK LVIEFSEEYP NKPPTVRFLS KMFHPNVYAD GSICLDILQN RWSPTYDVSS ILTSIQSLLD EPNPNSPANQ QAAQLYQENK REYEKRVSAI VEQSWNDS

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Ubc2/UBE2B protein under reducing conditions stained with Coomassie blue