

Human RUNX3/CBFA3 protein (Recombinant) (His-Tag) (STJP016716)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 0.1M NaCl, 1mM DTT
Molecular Weight	17.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	864
Gene Symbol	RUNX3
UniProt ID	RUNX3_HUMAN
Immunogen Region	53-186aa
Immunogen Sequence	MRSMDVDVLAD HAGELVRTDS PNFLCSVLPS HWRCNKTLPV AFKVVALGDV PDGTVVTVMA GNDENYS AEL RNASAVMKNQ VARFNDLRFV GRSGRGKSFT LTITVFTNPT QVATYHRAIK VTVDGPREPR RHRQK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RUNX3 CBFA3 protein under reducing conditions stained with Coomassie blue