

Human GRHPR protein (Recombinant) (His-Tag) (STJP016717)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.2M NaCl, 5mM DTT, 20% Glycerol
Molecular Weight	37.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9380
Gene Symbol	GRHPR
UniProt ID	GRHPR_HUMAN
Immunogen Region	1-328aa
Immunogen Sequence	MRPVRLMKVF VTRRIPAEGR VALARAADCE VEQWDSDEPI PAKELERGVA GAHGLLCLLS DHVDKRILDA AGANLKVIST MSVGIDHLAL DEIKKRGIRV GYTPDVLTD TAEAVSLLL TTCRRLPEAI EEVKNNGGWS WKPLWLCGYG LTQSTVGIIG LGRIGQAIAR RLKPFQVQRF LYTGRQPRPE EAAEFQAEFV STPELA

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human GRHPR protein 3µg under reducing conditions Coomassie blue stain