

Human RRM2 protein (Recombinant) (His-Tag) (STJP016722)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.1M NaCl
Molecular Weight	47 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6241
Gene Symbol	RRM2
UniProt ID	R1R2_HUMAN
Immunogen Region	1-389aa
Immunogen Sequence	MLSLRVPLAP ITDPQQQLQLS PLKGLSLVDK ENTTPALSGT RVLASKTARR IFQEPTEPKT KAAAPGVEDE PLLRENPRRF VIFPIEYHDI WQMYKAEAS FWTAEEDVLS KDIQHWESLK PEERYFISHV LAFFAASDGI VNENLVERFS QEVQITEARC FYGFQIAMEN IHSEMYSLI DTYIKDPKER EFLFNAIETM PCVKKK

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human RRM2 protein 3µg under reducing conditions stained with Coomassie blue