

Human NCK1 protein (Recombinant) (His-Tag) (STJP016727)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% Glycerol, 1mM EDTA, 50mM NaCl
Molecular Weight	45 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4690
Gene Symbol	NCK1
UniProt ID	NCK1_HUMAN
Immunogen Region	1-377aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MAEEVVVAK FDYVAQQEQE LDIKKNERLW LLDDSKSWWR VRNSMKNKTGF VPSNYVERKN SARKASIVKN LKDTLGIGKV KRKPSVPDSA SPADDSFVDP GERLYDLNMP AYVKFNYMAE REDELSLIKG TKVIVMEKCS DGWWRGSYNG QVGWFP SNYV TEEGDSPLGD HVGSLSEKLA AVVNNLNTGQ VLHVQALYP FSSSNDEE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human NCK1 protein 3µg under reducing conditions stained with Coomassie blue