

Human DCTD protein (Recombinant) (His-Tag) (STJP016755)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 1mM DTT, 20% Glycerol
Molecular Weight	22.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1635
Gene Symbol	DCTD
UniProt ID	DCTD_HUMAN
Immunogen Region	1-178aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MSEVSCCKRD DYLEWPEYFM AVAFLSAQRS KDPNSQVGAC IVNSENKIVG IGYNGMPNGC SDDVLPWRRRT AENKLDTKYP YVCHAEINAI MNKNSTDVKG CSMYVALFPC NECAKLIQA GIKEVIFMSD KYHDSDEATA ARLLFNMAV TFRKFIPKCS KIVIDFDSIN SRPSQKLQ

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



15% SDS-PAGE (3µg)

SDS-PAGE gel of recombinant human DCTD protein under reducing conditions Coomassie blue stain