

Human WDR5 protein (Recombinant) (His-Tag) (STJP016772)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 2mM DTT, 0.1M NaCl.
Molecular Weight	38.8 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11091
Gene Symbol	WDR5
UniProt ID	WDR5_HUMAN
Immunogen Region	1-334aa
Immunogen Sequence	MATEEKKPET EAARAQPTPS SSATQSKPTP VKPNYALKFT LAGHTKAVSS VKFSPNGEWL ASSSADKLIK IWGAYDGKFE KTISGHKLG I SDVAWSSDSN LLVSASDDKT LKIWDVSSGK CLKTLKGHSN YVFCCNFPQ SNLIVSGSFD ESVRIWDVKT GKCLKTLPAH SDPVSAVHFN RDGSLIVSSS YDGLCRIWDT ASGQCL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human WDR5 protein 3µg under reducing conditions Coomassie blue stain