

Human NUDT5 protein (Recombinant) (His-Tag) (STJP016773)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT, 0.1M NaCl
Molecular Weight	26.5 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11164
Gene Symbol	NUDT5
UniProt ID	NUDT5_HUMAN
Immunogen Region	1-219aa
Immunogen Sequence	MESQEPTSS QNGKQYIIE ELISEGKWKV LEKTTYMDPT GKTRTWESVK RTRRKEQTAD GVAVIPVLQR TLHYECIVLV KQFRPPMGY CIEFPAGLID DGETPEAAAL RELEEETGYK GDIAECSPAV CMDPGLSNCT IHIVTVTING DDAENARPKP KPGDGEFVEV ISLPKNDLLQ RLDALVAEEH LTVDARVYSY ALALKH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human NUDT5 protein under reducing conditions Coomassie blue stain