

Human Snurportin 1/SNUPN protein (Recombinant) (His-Tag) (STJP016780)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol 0.1M NaCl
Molecular Weight	43.3 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10073
Gene Symbol	SNUPN
UniProt ID	SPN1_HUMAN
Immunogen Region	1-360aa
Immunogen Sequence	MEELSQLALAS SFSVSQLDLS TAAPHRLSQ YKSKYSSLEQ SERRRRLLEL QKSKRLDYVN HARRLAEDDW TGMSEEEENK KDDEEMDIDT VKKLPKHIAN QLMLSEWLID VPSDLGQEWI VVVCVPGKRA LIVASRGSTS AYTKSGYCVN RFSSLLPGGN RRNSTAKDYT ILDCIYNEVN QTYVLDVVC WRGHPFYDCQ TDFRFY

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Snurportin 1/SNUPN protein under reducing conditions stained with Coomassie blue