

Human GALK1 protein (Recombinant) (His-Tag) (STJP016781)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol 0.1M NaCl, 1mM DTT
Molecular Weight	44.4 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2584
Gene Symbol	GALK1
UniProt ID	GALK1_HUMAN
Immunogen Region	1-392aa
Immunogen Sequence	MAALRQPQVA ELLAEARRAF REEFGAPEL AVSAPGRVNL IGEHTDYNQG LVLPMALELM TVLVGSPRKD GLVSLTTSE GADEPQRLQF PLPTAQRSLE PGTPRWANYV KGVIIQYYPAA PLPGFSAVVV SSVPLGGGLS SSASLEVATY TFLQQLCPDS GTIAARAQVC QQAHSFAGM PCGIMDQFIS LMGQKGHALL IDCRL

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human GALK1 protein under reducing conditions Coomassie blue stain