

E.coli nfsB protein (Recombinant) (His-Tag) (STJP016788)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 50mM NaCl, 1mM DTT, 10% Glycerol
Molecular Weight	26.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	1-217aa
Immunogen Sequence	MDIISVALKR HSTKAFDASK KLTPEQAEQI KTLQYSPSS TNSQPWHFIV ASTEEGKARV AKSAAGNYVF NERKMLDASH VVVFCAKTAM DDVWLKLVVD QEDADGRFAT PEAKAANDKG RKFFADMHRK DLHDDAEWMA KQVYLVGNF LLGVAALGLD AVPIEGFDA ILDAEFLKE KGYTSLVVVP VGHHSVEDFN ATLPKS

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant E. coli nfsB protein under reducing conditions stained with Coomassie blue