

Human LC3A/MAP1LC3A protein (Recombinant) (His-Tag) (STJP016796)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% Glycerol, 0.1M NaCl
Molecular Weight	15.1 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84557
Gene Symbol	MAP1LC3A
UniProt ID	MLP3A_HUMAN
Immunogen Region	1-120aa
Immunogen Sequence	MPSDRPFKQR RSFADRCKEV QQIRDQHPSK IPVIERKYK EKQLPVLDKT KFLVPDHSVNM SELVKIIRRR LQLNPTQAFF LLVNHSMVS VSTPIADIYE QEKDEDEGFLY MVYASQETFG

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel showing recombinant human LC3A MAP1LC3A protein 3µg under reducing conditions Coomassie blue stain