

Human Aminoacylase/ACY1 protein (Recombinant) (His-Tag) (STJP016799)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 10% Glycerol
Molecular Weight	48kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at-20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	95
Gene Symbol	ACY1
UniProt ID	ACY1_HUMAN
Immunogen Region	1-408aa
Immunogen Sequence	TSKGPEEEH PSVTLFRQYL RIRTVPKPD YGAAVAFFEE TARQLGLGCQ KVEVAPGYV TVLTWPGTNP TLSSILLNSH TDVVPVFKKH WSHDPFEAFK DSEGYYIARG AQDMKCVSIQ YLEAVRRLKV EGHFRPRTIH MTFVPDEEVG GHQGMELFVQ RPEFHALRAG FALDEGIANP TDAFTVFYSE RSPWWWVRVTS TGRPGH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Aminoacylase ACY1 protein under reducing conditions Coomassie blue stain