

## Human Fascin protein (Recombinant) (His-Tag) (STJP016802)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 20% Glycerol, 2mM DTT, 100mM NaCl
Molecular Weight	56.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6624</a>
Gene Symbol	<a href="#">FSCN1</a>
UniProt ID	<a href="#">FSCN1_HUMAN</a>
Immunogen Region	1-493aa
Immunogen Sequence	MTANGTAEAV QIQFGLINCG NKYLTAEEAFG FKVNASASSL KKKQIWTLEQ PPDEAGSAAV CLRSHLGRYL AADKDGNTVC EREVPGPDCR FLIVAHDDGR WSLQSEAHRR YFGGTEDRLS CFAQTVSPA E KWSVHIAMHP QVNIYSVTRK RYAHLSARPA DEIAVDRDVP WGVDSLITLA FQDQRYSVQT ADHRFLRH DG RLVARP

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE of recombinant human Fascin protein 3µg under reducing conditions Coomassie blue stain