

## Human TALDO1 protein (Recombinant) (His-Tag) (STJP016806)

### GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

### PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol
Molecular Weight	39.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

Gene ID	<a href="#">6888</a>
Gene Symbol	<a href="#">TALDO1</a>
UniProt ID	<a href="#">TALDO_HUMAN</a>
Immunogen Region	1-337aa
Immunogen Sequence	MGSSHHHHH SGLVPRGSH MSSSPVQRQR MESALDQLKQ FTTVVADTGD FHAIDEYKPKQ DATTNPSLIL AAAQMPAYQE LVEEAIAAYGR KLGGSQEDQI KNAIDKLFVL FGAEILKKIP GRVSTEV DAR LSFDKDAMVA RARRLIELYK EAGISKDRIL IKLSSTWEGI QAGKELEEQH GIHCNMTLLF SFAQAVACAE AGVTLISPFV GRILDWHVAN TDKKSYEP

### ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human TALDO1 protein 3µg under reducing conditions with Coomassie blue staining