

Human OGG1 protein (Recombinant) (His-Tag) (STJP016808)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in phosphate-Buffered Saline (pH 7.4) containing 40% Glycerol
Molecular Weight	41.2 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4968
Gene Symbol	OGG1
UniProt ID	OGG1_HUMAN
Immunogen Region	1-345aa
Immunogen Sequence	MPARALL PRRMGHRTLA STPALWASIP CPRSELRLDL VLPSSGQSFWRV REQSPAHWSG VLADQVWTLT QTEQLHCTV YRGDKSQASR PTPDELEAVR KYFQLDVTLA QLYHHWGSVD SHFQEVAKF QGVRLLRQDP IECLFSFICS SNNNIARITG MVERLCQAFG PRLIQLDDVT YHGFPQLQAL AGPEVEAHLR KLGLGY

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human OGG1 protein 3µg under reducing conditions stained with Coomassie blue