

Human Annexin A8/ANXA8 protein (Recombinant) (His-Tag) (STJP016822)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES

Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 1mM DTT, 1mM EDTA, 10% Glycerol
Molecular Weight	39.0 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	653145/728113
Gene Symbol	ANXA8
UniProt ID	ANXA8_HUMAN
Immunogen Region	1-327aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MAWWKSWIEQ EGVTVKSSSH FNPDPDAETL YKAMKGIGTN EQAIIDVLTK RSNTQRQQIA KSFKAQFGKD LTETLKSELS GKFERLIVAL MYPPYRYEAK ELHDAMKGLG TKEGVIIIL ASRTKNQLRE IMKAYEEDYG SSLEEDIQAD TSGYLERILV CLLQGSRDDV SSFVDPGLAL QDAQDLYAAG EKIRGTDEMK FITILCTR

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human Annexin A8 ANXA8 protein 3 μ g under reducing conditions with Coomassie blue stain