

Human Galectin-8/LGALS8 protein (Recombinant) (His-Tag) (STJP016828)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE/Bioactivity
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 0.1M NaCl, 1mM DTT, 10% Glycerol
Molecular Weight	37.9 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3964
Gene Symbol	LGALS8
UniProt ID	LEG8_HUMAN
Immunogen Region	1-317aa
Immunogen Sequence	MMLSLNLLQN IYNPVIPIV GTIPDQLDPG TLIVIRGHVP SDADRFQVDL QNGSSMKPRA DVAFHFNPRF KRAQCIVCNT LINEKWGREE ITYDTPFKRE KSFEIVIMVL KDKFQVAVNG KHTLLYGHRI GPEKIDTLGI YGKVNHSIG FSFSSDLQST QASSLELLEI SRENVPKSGT PQLRLPFAAR LNTPMGPGRT VVVKGE

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human Galectin-8 LGALS8 protein under reducing conditions stained with Coomassie blue