

firefly Luciferase protein (Recombinant) (His-Tag) (STJP016840)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% Glycerol
Molecular Weight	63 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Immunogen Region	1-550aa
Immunogen Sequence	MEDAKNIKK GPAPFYPLED GTAGEQLHKA MKRYALVPGT IAFTDAHIEV NITYAEYFEM SVRLAEAMKR YGLNTNHRIV VCSENSLQFF MPVLGALFIG VAVAPANDIY NERELLNSMN ISQPTVVVFS KKGLQKILNV QKKLP IIQKI IIMDSKTDYQ GFQSMYTFVT SHLPPGFNEY DFVPESFDRD KTIALIMNSS GSTGLP

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant firefly luciferase protein under reducing conditions Coomassie blue stain