


Human ECHS1 protein (Recombinant) (His-Tag) (STJP016843)

GENERAL INFORMATION	
Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli
PRODUCT PROPERTIES	
Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 20% Glycerol, 100mM NaCl
Molecular Weight	30.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	1892
Gene Symbol	ECHS1
UniProt ID	ECHM_HUMAN
Immunogen Region	28-290aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MASGANFEYI IAEKRGKNNT VGLIQLNRPK ALNALCDGLI DELNQALKIF EEDPAVGAIV LTGGDKAF AA GADIKEMQNL SFQDCYSSKF LKHWDH LTQV KKPVIAAVNG YAFGGGCELA MMCDI IYAGE KAQFAQPEIL IGTIPGAGGT QRLTRAVGKS LAMEMVLTGD RISAQDAKQA GLVSKICPVE TLVEEAIQCA EKIASNSK
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human ECHS1 protein 3µg under reducing conditions with Coomassie blue stain