

Human GMPR2 protein (Recombinant) (His-Tag) (STJP016844)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 1mM DTT, 10% Glycerol
Molecular Weight	40 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51292
Gene Symbol	GMPR2
UniProt ID	GMPR2_HUMAN
Immunogen Region	1-348aa
Immunogen Sequence	MPHIDNDVKL DFKDVLLRPK RSTLKSRSSEV DLTRSFSEFRN SKQTYSGVPI IAAANMDTVGT FEMAKVLCKF SLFTAVHKHY SLVQWQEFAG QNPDCLEHLA ASSGTGSSDF EQLEQILEAI PQVKYICLDV ANGYSEHFVE FVKDVRKRFP QHTIMAGNVV TGEMVEELIL SGADIIVGI GPGSVCTTRK KTGVGYPQLS AVMECA

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel showing recombinant human GMPR2 protein 3µg under reducing conditions with Coomassie blue stain