

Human ANP32A protein (Recombinant) (His-Tag) (STJP016852)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT, 100mM NaCl, 0.1mM PMSF
Molecular Weight	30.7 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8125
Gene Symbol	ANP32A
UniProt ID	AN32A_HUMAN
Immunogen Region	1-249aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MEMGRRIHLE LRNRTPSDVK ELVLDNSRSN EGKLEGLTDE FEELEFLSTI NVGLTSIANL PKLNKLLKLE LSDNRVSGGL EVLAEKCPNL THLNLSGNKI KDLSTIEPLK KLENLKSIDL FNCEVTNLND YRENVFKLLP QLTYLDGYDR DDKEAPDSA EGYVEGLDDE EDEDEEEYD EDAQVVEDEE DEDEEEGEE EDVSGEEE

ADDITIONAL INFORMATION

Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.
------	---



SDS-PAGE gel of recombinant human ANP32A protein under reducing conditions