

Human SEPT2/SEPTIN2 protein (Recombinant) (His-Tag) (STJP016857)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.25 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 40% Glycerol, 0.1M NaCl, 1mM DTT
Molecular Weight	43.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4735
Gene Symbol	SEPTIN2
UniProt ID	SEPT2_HUMAN
Immunogen Region	1-361aa
Immunogen Sequence	MSKQQPTQFI NPETPGYVGF ANLPNQVHRK SVKKGFEFTL MVLGSEGLGK STLINSLFLT DLYPERVISG AAEKIERTVQ IEASTVEIEE RGVKLRLLTVV DTPGYGDAIN CRDCFKTIIS YIDEQFERYL HDESGLNRRH IIDNRVHCCF YFISPFHGHL KPLDVAFMKA IHNKVNIVPV IAKADTLTLK ERERLKKRIL DEIEEH

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human SEPT2/SEPTIN2 protein under reducing conditions stained with Coomassie blue