


Human DUTPase/DUT protein (Recombinant) (His-Tag) (STJP016865)

GENERAL INFORMATION	
Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli
PRODUCT PROPERTIES	
Concentration	1 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 8.0) containing 10% Glycerol, 1mM DTT, 0.1M NaCl.
Molecular Weight	21.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.
TARGET INFORMATION	
Gene ID	1854
Gene Symbol	DUT
UniProt ID	DUT_HUMAN
Immunogen Region	70-252aa
Immunogen Sequence	MGSSHHHHHH SSGLVPRGSH MASTVGAAGW KGELPKAGGS PAPGPETPAI SPSKRARPAE VGGMQLRFAR LSEHATAPTR GSARAAGYDL YSAYDYTIPT MEKAVVKTDI QIALPSGCGY RVAPRSGGLAA KHFIDVGAGV IDEDYRGNVG VVLFNFGKEK FEVKKGDRIA QLICERIFYP EIEEVQALDD TERGSGGFGS TGKN
ADDITIONAL INFORMATION	
Note	STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human DUTPase DUT protein under reducing conditions stained with Coomassie blue