

Human DCK protein (Recombinant) (His-Tag) (STJP016878)

GENERAL INFORMATION

Product Type	Proteins
Applications	SDS-PAGE
Host / Source	E.coli

PRODUCT PROPERTIES


Concentration	0.5 mg/mL
Formulation	Liquid in 20mM Tris-HCl buffer (pH 7.5) containing 1mM DTT, 0.1mM PMSF, 2mM EDTA, 10% Glycerol
Molecular Weight	34.6 kDa
Storage Instruction	For short term storage, keep at +2C to +8C for up to 1 week. For long term storage, aliquot and store at -20C, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1633
Gene Symbol	DCK
UniProt ID	DCK_HUMAN
Immunogen Region	1-260aa
Immunogen Sequence	MATP PKRSCPSFSA SSEGTRIKKI SIEGNIAAGK STFVNILKQL CEDWEVVPEP VARWCNVQST QDEFEELTMS QKNGGNVLQM MYEKPERWSF TFQTYACLSR IRAQLASLNG KDKAEKPV LFFERSVYS DR YIFASNLYES ECMNETEW TI YQDWHDW MN QFGQSLELDG IYQLQATPET CLHRIY

ADDITIONAL INFORMATION

Note STRICTLY FOR FURTHER SCIENTIFIC RESEARCH USE ONLY (RUO). MUST NOT TO BE USED IN DIAGNOSTIC OR THERAPEUTIC APPLICATIONS.



SDS-PAGE gel of recombinant human DCK protein 3µg under reducing conditions stained with Coomassie blue